





InterSpec gamma radiation analysis software

William Johnson 201810XX


wcjohns@sandia.gov

<https://github.com/sandialabs/InterSpec>



- InterSpec helps to interactively perform nuclide identification using gamma data from a wide variety of detectors
 - Uses a peak-based analysis approach with comprehensive nuclide and shielding databases
 - Can be used with identiFINDER, RadSeeker, other NaI, HPGe, LaBr, and more detectors
- Can be used to perform nuclide activity quantification, determine amount and/or type of shielding, nuclide age
 - What information can be extracted depends on nuclides present, shielding amounts, and data quality

4 InterSpec Analysis Overview



Analyzing data in InterSpec primarily happens by fitting for “photo-peaks” in gamma spectrum



Nuclide: Show Lines **Persist**

Age: $\lambda=1.2e+03$ y **Clear All**

Min Amp:

Detector: NaI 3x3

Thickness

Gammas
 X Rays
 Alphas
 Betas

Energy (keV)	B.R.	Parent	Mode	Particle
47.7340	4.97e-05	Ho166m	xray	xray
48.2210	0.108	Ho166m	xray	xray
49.1280	0.192	Ho166m	xray	xray
55.4800	0.0201	Ho166m	xray	xray
55.6740	0.0387	Ho166m	xray	xray
56.0540	0.00074	Ho166m	xray	xray
57.1420	0.013	Ho166m	xray	xray
57.3130	0.00385	Ho166m	xray	xray
73.4500	0.0001452	Ho166m	β^-	γ

Since the expected yields of gammas at a given energy is known, as well as the efficiency of detection for a given detector at each energy is known – InterSpec can use this information to back-out: source strength, shielding amount, possibly shielding type, possibly nuclide age, and more information

To use InterSpec to analyze spectra, you will generally perform the following steps:


- Make sure the spectrum has a good energy-calibration. This usually done using a background, or a known source, taken at around the same time as the spectrum of interest
- Identify the source nuclides (or x-rays, or nuclear reactions) that cause each of the peaks in the spectrum
- Fit for the nuclides activities, maybe ages, and shielding amount and types

Energy Calibration Check: Spectrum of a known Cs137 source



Expected Energy
(see "Reference Photopeaks" tab)

Observed Peak Energy



Usually, you check that one or a few peaks, that you know what energy they should be, are actually at those energies.

Energy Calibration Check: fixing calibration up




Since we know what energies the peaks should be at, we can have InterSpec fit for the correct energy calibration coefficients

The screenshot shows the InterSpec software interface. On the left, there are input fields for fit terms: Offset Term (0.000000), Linear Term (3.674617), and Quad. Term (0.000000). Below these are checkboxes for 'fit' and 'Apply to: Foreground, 2nd Spec, Background'. On the right, a table displays calibration data for Cs137. A red circle highlights the table header and the first row. A black arrow points from the text on the left to the 'Fit Coeffs' button at the bottom right, which is also circled in red.

Nuclide	Calib. Peak	Mean	Photopack	Difference
Cs137	<input checked="" type="checkbox"/> true	791.59	661.66 keV	-129.93 keV

Or you can graphically recalibrate by right-click dragging the peaks to where they are supposed to be



Or you can manually play with the calibration coefficients

Offset Term $\times 10^0$ keV fit

Linear Term $\times 10^3$ keV/ F_{WF} fit

Quad. Term $\times 10^0$ keV/ F_{WF}^2 fit


Apply to: Foreground 2nd Spec Background

9 Nuclide Identification:



InterSpec uses the photo-peaks in the spectrum to help determine which nuclides are present

Photopeak's observed in the spectrum




The expected photopeak locations for Ba133

Nuclide Identification (cont):



InterSpec does not perform automatic nuclide identification, but provides tools to help you determine them




Clicking on a trial nuclide shows the photopeak reference lines for that nuclide – this allows determining which nuclide most closely matches observed data

Here we are searching for nuclides that have the three photopeaks highlighted

Nuclide Identification (cont):



Shielding or nuclide age can alter what photopeak are observed in the data – InterSpec can help with this



Spectrum Files **Peak Manager** **Reference Photopeaks** **Energy Calibration** **Nuclide Search**

Nuclide: Show Lines

Age: $\lambda=1.6e+03$ y

Min Amp: Prompt

Detector: IdentIFINDER-NGH Gammas X Rays

Alphas Betas

 $\rho=11$ g/cm³ Thickness

Energy (keV)	B.R.
10.1370	0.000156
10.1720	1.496e-11
10.2680	1.337e-10
10.4500	1.654e-06
10.5510	1.464e-05
10.7310	0.009419
10.8390	0.08407
10.9940	1.921e-12
11.0160	0.0003819


Entering shielding amount, nuclide age, or detector response changes amplitude of reference photopeaks

(calculating shielding/age from data to be shown later in presentation)


Nuclide and Shielding Quantification:



- InterSpec uses the amplitude of peaks, in combination with the known emission rate of gamma radiation, a known distance, and the detector response function (efficiency of detecting gamma at a given energy) to calculate source strength
- If an unknown amount of shielding is present: the amount of shielding, and possibly its effective atomic number, may be able to be calculated from the data
 - The probability of gammas interacting with the shielding is energy and atomic number dependent; InterSpec will use the relative amplitudes of peaks at different energies in the spectrum to calculate how much and what type of shielding there is




InterSpec allows calculating activities of multiple nuclides



If multiple nuclides contribute to a peak, InterSpec will account for this

Nuclide Quantification: Fitting for values



The “Activity/Shielding Fit” tool uses the peaks you fit for to then fit for source strength/age and shielding

Using InterSpec

InterSpec requirements:



- InterSpec works on all major platforms
 - Using on iPhone/iPad and Android works surprisingly well!
- Open-source: <https://github.com/sandialabs/InterSpec>
 - You can make any changes you want to it, and/or compile it from source
- Does not require installing on Windows – just unzip the package and run the EXE
 - On the google play app store, and we are (slowly) working on adding to Apple app stores – currently distributed for iOS ad-hoc
- InterSpec accepts data from nearly all commercially available RIID detectors and detection systems
- InterSpec comes with a few generic detector response functions, but ones for specific model detectors will have to be given you separately






We will now perform a few example analysis – the next 23 slides is are meant for later reference

Example 1




A metal container was located that was giving off increased radiation

- A 5 minute spectrum was taken from 1 m away, along with a 5 minute representative background; both using the same 3x3 NaI detector†
- You would like to determine:
 - What nuclides are present inside the box?
 - The nuclide(s) activity?
 - How much shielding is present?



† 7.62 cm diameter by 7.62 cm long NaI crystal - will produce the conceptually same spectroscopic data as an identiFINDER or RadSeeker


Example 1 (cont):



Easiest way to load spectra is to drag-n-drop from the operating system


- You can display up to three spectra at a time, but peak-fitting is limited to the foreground
- If you have previously worked with the spectra, you may be prompted if you want to pick up where you left off (you can also explicitly save state to InterSpecs database, or “tag” or “branch”)

Example I (cont):



- The background peaks at 1460 keV (K40), 609 keV (Ra226), and 2614 keV (Th232) look to be at the expected energies
- The energy calibration is acceptable, and it looks like we got the correct files off the detector

Example I (cont):




There are two obvious photo-peaks in your item of interest spectrum that are not in the background spectrum.

You can fit peaks by double-clicking near them.

Example I (cont):


Lets figure out what nuclide causes these photopeaks by using the “Nuclide Search” tab

When on the “Nuclide Search” tab if you click on a place in the spectrum – that energy is filled in to search on



The number of energies can be changed by using the + and - icons


Example I (cont):



Clicking on a nuclide will cause the reference photopeak lines for that nuclide to be shown

- Here Fm256 is clearly not right because we would see many other peaks

Example I (cont):




Its usually best to sort this table by the “Branch Ratio” column (default), but other ordering can be changed by clicking on the headers for other columns

It’s a little hard to see, but the reference lines for Na22 line up exactly with the observed peaks – and we aren’t missing any peaks we would expect for Na22 – so this is our nuclide

Example I (cont):

Now we need to assign nuclides to these peaks:



If you right-click on the peak, and go to the “Change Nuclide” menu-item, you can select Na22, or there are a number of other possible isotopes that it could be from

Example I (cont):

You can also assign nuclides to peaks, in a number of other ways:

By typing the nuclide name into the “Peak Manager” tab

Spectrum Files		Peak Manager	Refer
Nuclide	Mean	FWHM	
Na22	509.26	35.61	
Na22	1272.51	60.02	

The “Peak Editor” tool (opened by right clicking on peak and selecting “Peak Editor...”)

Parameter	Value	Uncertainty	Fit
Centroid	509.2592	0.1805	<input checked="" type="checkbox"/>
FWHM	35.6116	0.4598	<input checked="" type="checkbox"/>
Amplitude	20699.3008	217.3821	<input checked="" type="checkbox"/>
χ^2/DOF	4.7400	0.0000	
ROI Start (keV)	433.1371	0.0000	
ROI End (keV)	563.4930	0.0000	
Cont. P0	365.2334	3.0797	<input checked="" type="checkbox"/>
Cont. P1	-1.7892	0.0340	<input checked="" type="checkbox"/>
Nuclide	Na22		
Photopeak Label	510.9989 keV $I=1.0e+C$		
Peak Type	Gaussian		
Continuum	Linear		
Skew Type	None		

Reference Photopeaks

Nuclide: Show Lines

Age: $\lambda=2.60$ y


Min Amp:

Detector: NaI 3x

Gammas X Rays Alphas Betas

Thickness:

Or if you are showing reference photo-peak lines when you fit for a peak, the peak will automatically be associated with that nuclide





Mousing over a peak causes a pop-up that shows info about the peak, including nuclide


Example I (cont):



Fit for the source activity using the “Activity/Shielding Fit” tool



Example I (cont):



Select the detector response function you want to use by clicking here



Activity/Shielding Fit

(observed-model)/uncert

Energy (keV)

<dev>=0.089σ

Distance: 100 cm

Add Shielding: Material Generic

Fe

$\rho=7.9 \text{ g/cm}^3$, $\overline{AN}\sim 26.0$

Thickness 3.40 (± 0.23) cm Fit

Nuclide	Mean	Photo...	Use
Na22	509.99	511.00 keV	<input checked="" type="checkbox"/> true
Na22	1272.66	1274.54 keV	<input checked="" type="checkbox"/> true

Nuclide	Activi...	Fit Act.	Age	Fit Age
Na22	39.42 uCi	<input checked="" type="checkbox"/> true	NA	

All Peaks

Multiple nuclides contribute to peak

Subtract Background Peaks

Isotopes of same element same age

Show Chi Graphic

Perform Model Fit

Close

Enter distance of measurement

Example I (cont):



Select which peaks you would like to use

Activity/Shielding Fit

Detector: NaI 3x3

Distance: 100 cm

Add Shielding: Material Generic

Fe
 $\rho=7.9 \text{ g/cm}^3$, $\overline{AN}\sim 26.0$
Thickness 3.40 (± 0.23) cm Fit

Perform Model Fit

Close

Energy (keV) $\langle dev \rangle = 0.089\sigma$

(observed-model)/uncert


Nuclide	Mean	Photo...	Use
Na22	509.99	511.00 keV	<input checked="" type="checkbox"/> true
Na22	1272.66	1274.54 keV	<input checked="" type="checkbox"/> true

All Peaks

Multiple nuclides contribute to peak
 Subtract Background Peaks
 Isotopes of same element same age
 Show Chi Graphic

Nuclide	Activi...	Fit Act.	Age	Fit Age
Na22	39.42 uCi	<input checked="" type="checkbox"/> true	NA	

Example I (cont):



Add a shielding

Here we added a
single iron shielding

Make sure you are fitting for
activity, and shielding thickness



Activity/Shielding Fit

(observed-model)/uncert

Energy (keV)

<dev>=0.089 σ

Detector: NaI 3x3

Distance: 100 cm

Add Shielding: Material Generic

Fe

$\rho=7.9 \text{ g/cm}^3$, $\overline{AN}\sim 26.0$

Thickness 3.40 (± 0.23) cm Fit

Nuclide	Mean	Photo...	Use
Na22	509.99	511.00 keV	<input checked="" type="checkbox"/> true
Na22	1272.66	1274.54 keV	<input checked="" type="checkbox"/> true

Nuclide	Activi...	Fit Act.	Age	Fit Age
Na22	39.42 uCi	<input checked="" type="checkbox"/> true	NA	

All Peaks

Multiple nuclides contribute to peak

Subtract Background Peaks

Isotopes of same element same age


Show Chi Graphic

Perform Model Fit

Close


Click here to fit for desired quantities

Example I (cont)



We fit for 40 uCi of Na22, with 3.4 cm of Iron Shielding


- Truth level: 50 uCi at 3.48 cm Fe - not bad!



This “Chi2” graphic shows you how many statistical sigmas the fit peak areas are off from what is predicted for the fit activity and shielding


- Here we fit for two quantities (activity and shielding thickness), using two peaks, so both peak areas should be at essentially zero sigmas. If we had more peaks we would expect a distribution ranging between about ± 5 sigmas (or maybe a little more)

Example 2:



For this problems we are told the same 3x NaI detector was used, but no background was provided, and shielding is lead

Example 2 (cont):




Truth is: 45 uCi Eu152, 25 uCi Cs137, with 0.1765cm Pb Shielding

Example 3:



A chunk of Thorium metal is measured at 1m, using a 40% HPGe detector – determine activity.



Example 3 (cont): Truth: 2.54 cm sphere of Thorium



Some things to note:

- Thorium self attenuates – InterSpec can account for this, but always assumes source is a sphere
- X-rays present in spectrum implying little or no shielding – InterSpec also fit for no additional shielding outside Th
- Single and double escape peaks, and 511 keV peak.

InterSpec isn't magic:



- If you have a source with a single peak, you can only fit for a single quantity: Activity, or Shielding amount. You can't fit for shielding type (atomic number), or nuclide age.
- Fitting for shielding type, it is best to have multiple peaks over a wide energy region. Also, fitting for shielding type is best done for HPGe detectors, with sources that have many peaks.
- If you can fit for a nuclide age depends on if the gamma signature of the nuclide changes over time
- It assumes a qualified analyst is performing the analysis. Lots of things can be messed up, like assigning the wrong nuclide to a peak, not realizing a nuclide needs to be aged, etc



- It is currently a “in my personal time” project, but we are looking for funding sources
- It does have a few hundred unit tests, as well as a end-to-end testing mechanism (none of these are in the publicly available code), but it is, in general, not strictly validated. The assumption is that the qualified analyst using the program could catch any issues/errors.
 - But if you do find any bugs, let me know – I would love to fix them!
- The help documentation needs a lot more work

Other useful things in InterSpec



Spectrum

 Foreground
 Secondary
 Background

Allow Edit

 Yes
 No

File Information

File Name:	ba133_source_640s_20100317.n42	Mem Size:	581.0 kb
Inspection:	Lane: 1	Location: Example Location	RIID Analysis
Instru. Type: SpecPortal	Manufacturer: ORTEC	Model: OSASP	
Instru. ID: Serial #1XXX	UUID: e4592733-68e0-4f07-b1c3-0355739c283c		

File Remarks:
 Measured Data
 Occupancy number = 18
 Local End Time 2010-04-17T11:43:59.159
 Local Start Time 2010-04-17T11:33:19.409
 DNDORadiationMeasurement

Measurement Information



Date/Time: 2010-Apr-17 17:33:19.409000	Live Time: 629.17 s	Real Time: 639.92 s	
Det. Name: A1	Sample Num: 1	Energy(keV): -0.2 to 3084.4	Num Channels: 16384
Sum Gamma: 305570.00	Gamma CPS: 485.672	Sum Neutron: N/A	Neutron CPS: N/A
Latitude: dec or deg min' sec" N/S	Longitude: dec or deg min' sec" E/W	Position Time:	<input type="button" value="Show Map"/>

Description: **Source Type:**

Spectra Remarks:

◀ 1 of 8 ▶

Spectrum files often contain lots of extra information that are useful to figure things out. InterSpec can show you a lot of this information, including: GPS coordinates, date/time, RIID analysis results, user entered notes, serial numbers, and more



Flexible nuclide decay information (including export to CSV), reference decay information, and decay calculation

- The nuclide database is quite comprehensive



Dose

Activity

Distance

Shielding

Inputs

Source: Gamma Neutron

Nuclide: Cs137 λ=30.07 y

Age: NA

Activity: 100 μCi

Distance: 100 cm

Shielding:
 Thickness 1.0 cm

Answer

units: rem/hr

29.60 urem/hr

Stay Time	
1.5 mrem <small>Dental X-Ray</small>	50.67 h
620 mrem <small>Typ. Yearly Background</small>	2.39 y
5 rem <small>Annual Occ. Limit</small>	19.27 y

Close

Dose

Activity

Distance

Shielding

Inputs

Source: Gamma Neutron

Nuclide: Cs137 λ=30.07 y

Age: NA

Dose: 100 μR/hr

Distance: 100 cm

Shielding:
 Thickness 1.0 cm

Answer

units: curries

337.81 uCi

Stay Time	
1.5 mrem <small>Dental X-Ray</small>	15.00 h
620 mrem <small>Typ. Yearly Background</small>	258.33 d
5 rem <small>Annual Occ. Limit</small>	5.70 y

Close

- Either calculate expected dose from a nuclide with a given activity, shielding, and distance
- Or if you have the dose (like from a pager-style detector, or identiFINDER, etc), you can calculate activity
- And similar for distance/shielding