



Sandia
National
Laboratories

InterSpec Familiarization

Will Johnson – wcjohns@sandia.gov

20220420

Help/Support: InterSpec@sandia.gov




Sandia National Laboratories is a multimission laboratory managed and operated by National Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.

SAND2022-4588 TR

Today's goals:

- Get you comfortable with opening and viewing spectrum files in InterSpec
- Get you familiar with basic peak-fitting, nuclide ID, energy calibration
- Have you be able to fit activity and shielding for a spectrum



Please work-along, or tinker around with InterSpec as we go along today

The example spectrum files at <https://sandialabs.github.io/InterSpec/tutorials/>



Some general information



This presentation uses InterSpec v1.0.10_rc3, but v1.0.9 is acceptable

- Windows, Linux, macOS: <http://github.com/sandialabs/InterSpec/releases/>
 - Windows Install: unzip downloaded file, and place resulting directory where ever you would like (e.x., C:\Program Files\InterSpec); manually pin to start-menu and associate file-types
- iPhone, iPad, Android, macOS: search the app-store for "InterSpec"

InterSpec is completely open-source (LGPL v2.1), and code available from:

- <https://github.com/sandialabs/InterSpec>
- <https://github.com/sandialabs/SandiaDecay>
- <https://github.com/sandialabs/SpecUtils>

Support, bug reports, and requests: InterSpec@sandia.gov

Today's presentation should take ~1 hour




Introduction to InterSpec


general use



InterSpec Capabilities:



Quickly view spectrum files and fit peaks





Perform nuclide ID



Determine nuclide activities, shielding amounts, isotopics, and ages

InterSpec Capabilities (continued):



Perform nuclide decay calculations, get reference information, export inf as CSV, etc

1/r² Calculator window. Use two measurement at different locations to find distance to an unseen source. E.g. when the source is behind a wall.

Near Measurement Intensity: 100

Far Measurement Intensity: 13

Background Intensity (optional): 1


Distance between measurements: 2

Power Law: Low Scatter or using Peak Area, 1/r²

Dist. near measurement to source: 1.068

Use the same units for near, background, and far measurements. Results are in same units used for distance between measurements.

1/r² calculator



Develop detector response functions (efficiency + FWHM)

InterSpec Capabilities (continued):



Stay Time	
1.5 mrem	66.94 h
620 mrem	3.16 y
5 rem	25.45 y

Flexible gamma dose calculations
 Activity \leftrightarrow Dose
 Activity \leftrightarrow Distance
 Etc.

Energy (keV)	Peak CPS	Flux ($\gamma/cm^2/s$)	$\gamma/4n/s$
122.17	1.363 ± 0.1	0.001461 ± 0.0	1.836e+04 ± 1531.7
244.83	1.377 ± 0.1	0.001721 ± 0.0	2.162e+04 ± 1693.2
344.04	6.401 ± 0.1	0.01078 ± 0.0	1.355e+05 ± 2792.7
411.81	0.5318 ± 0.0	0.001086 ± 0.0	1.365e+04 ± 23.7
443.22	0.5317 ± 0.0	0.001179 ± 0.0	1.482e+04 ± 22.4
778.98	2.02 ± 0.1	0.008339 ± 0.0	1.048e+05 ± 4058.0
866.11	0.5508 ± 0.1	0.002544 ± 0.0	3.197e+04 ± 2960.2
964.39	1.845 ± 0.1	0.009523 ± 0.0	1.197e+05 ± 4315.1
1086.17	1.512 ± 0.2	0.0088 ± 0.0	1.106e+05 ± 15948.6
1112.41	1.763 ± 0.2	0.01051 ± 0.0	1.321e+05 ± 15543.2
1409.47	2.336 ± 0.1	0.01265 ± 0.0	2.218e+05 ± 5284.8

Gamma flux calculator

Convert between radiation related units.
 Ex: 5 MBq, 2 nCi, 1.2rad, 15E-3gy, 0.2mrem, 8feet, 9milli-sievert

Input: 5 MBq
 Output: 135.1 uCi

Convert radiation relevant units

Total att. cross section	0.1434	cm²/g
Compton	0.1427	cm²/g
Rayleigh	0.002313	cm²/g
Photoelectric	0.000608	cm²/g
Pair production	0	cm²/g
Mass avrg atomic num	6.58	
Attenuation (optional):		
Density:	0.65	g/cm³
Thickness:	1 cm	
Trans. Frac.	0.911	
Intrinsic Efficiency	0.7643	
Solid Angle Fraction	0.1442	
Detection Efficiency	0.1102	
Total Efficiency	0.1004	

Cross-section calculator

And more!



(and also more under active development)

What InterSpec is not:




InterSpec is not a magic easy button – all features are interactive – doesn't tell you the answer automatically, but assists the user to perform analysis

Not full-spectrum analysis



Does not do data acquisition - only accepts data already taken

Does not generate reports or adhere to standardized methodologies



Initial analysis example (about 4 minutes)



InterSpec View Tools Help

Counts

0 100 200 300 400 500

2500 2600 Energy (keV)

Spectrum Files Peak Manager

Nuclide Mean FWHM

License and Terms

More in depth information

show at start when no spectra

Close

Search for Peaks

Peak: Add...

CSV

Welcome To InterSpec

You can use the references on the left to become more familiar with InterSpec, or you can pick up from a previous session or spectrum below. Alternatively, you can also drag and drop your own spectrum file onto InterSpec.

Saved States Example Spectra

- Ba-133 (16k bin N42)
- Passthrough (16k bin ICD1, 8 det., 133 samples)
- Background (16k bin N42)

10



Result: 46 uCi of Ir-192 with 6.20 mm of Pb shielding

Truth: 51 uCi of Ir-192 with 6.35 mm of Pb shielding






How to use InterSpec



Interacting with spectra



Loading Spectra




To load a spectrum file: drag-n-drop from your computer

- When file is over InterSpec, you can choose to open it as Foreground, Background, or Secondary

The **InterSpec**→**Open File...** menu option lets you browse for files



Loading Spectra (continued)




You can double-click on spectrum files in the operating system to open them in InterSpec

Window: you have to manually associate file types with InterSpec first

macOS/iOS/Android: spectrum files are automatically associated with InterSpec

Interacting with spectra



Zoom in: click and drag to the right - let go of the mouse to zoom in.


Zoom out: click and drag to left; zooming out happens continuously

To pan left or right: click and drag with right mouse-button. Or click and drag x-axis labels

Mouse wheel: up/down zooms in; left/right pans energy

Touch screens: one finger pans left/right, two horizontal fingers zoom in/out


Interacting with spectra (cont)



View → "Show Energy Slider" will bring up a strip chart



Basic peak-fitting



To fit a peak, double click near it.

Double clicking again will add another peak.


You can adjust the Region Of Interest (ROI), by moving mouse to edge of ROI

Right-click on a peak to get a menu to add or remove a peak from a ROI

Holding control and dragging across a region lets you fit for multiple peaks in a ROI




Peak Fitting (cont)



There is also an automated peak search that works decent, especially for HPGe

- Gets you to the ~80% level for HPGe spectra – e.g., may require some “clean-up”
- Does not use the currently showing reference lines, or detector response function (DRF) to identify peaks

Peak Fitting (cont)



You can also right-click on a peak and choose “Peak Editor...” from the popup menu to have more control over fitting a peak.

You can fix some of the peak quantities, and fit for other ones, manually adjust some of the parameters, select a continuum type, and more




If you forget how to do something:


Help → Welcome → Controls

or

Help → "Help Contents..."

or


There is often a  icon in lower-left of tools you can click on



The screenshot displays the InterSpec software interface. A help window is open, showing the 'Zoom-In/Out' tool. The window contains a small plot of a spectrum with a peak highlighted. Below the plot, the text reads: 'Left-Click Drag right to zoom-in. Left-Click Drag left to zoom-out.' The help window also includes a 'Welcome' button, a 'Controls' button, and a 'Close' button. The main interface shows a spectrum plot with 'Counts' on the y-axis and 'Energy (keV)' on the x-axis. The plot shows a peak at approximately 400 keV. The help window is titled 'Zoom-In/Out' and includes a list of tools: 'Zoom-In/Out', 'Pan Energy', 'Energy Slider', 'Peak Fit', 'ROI Fit', 'Peak Editor', 'Energy Calibration', 'Delete Peak', 'Adjust ROI', 'Scale Y', 'Count Gammas', and 'Keyboard Shortcuts'. The main interface also shows a 'Spectrum Files' tab, a 'Peak Manager' tab, and a 'Reference' tab. The bottom of the interface has a search bar and a 'Peak: Add...' button.



Important concept!




Spectrum Files		Peak Manager		Reference Photopeaks	Energy Calibration	Nuclide Search	
Nuclide	Mean	FWHM	Area	CPS	Photopeak	Differ...	
Ba133	79.54	1.16	22002	4.378 \pm 0.074	79.61 keV	0.07 keV	
Ba133	80.92	1.18	277536	55.22 \pm 0.12	81.00 keV	0.07 keV	
Ba133	276.22	1.42	41975	8.352 \pm 0.045	276.40 keV	0.18 keV	

Peaks can get associated with a nuclide, x-ray, or reaction

- The peak is actually associated with a specific gamma-line of the nuclide/x-ray/reaction



Associating peaks with a nuclide/reaction/x-ray



Associating a peak with a nuclide/reaction/x-ray:

- Keep track of nuclide ID
- Enables easy energy calibration, nuclide activity, isotope age, shielding, detector response function, calculate fluxes, ...


Associating a peak with a nuclide/x-ray/reaction – best method




Show reference lines **before** fitting for peaks

- Association will automatically be made when you fit the peak if there is a gamma-line near peak
- Peak color will be the same as the reference line, to help track ID

Associating a peak with a nuclide/x-ray/reaction – other methods





Right-click peak and going to “Change Nuclide” item will bring up some suggested nuclides




Tables listing peak information are all directly editable

The peak-editor also lets you select the nuclide and energy





By default peaks are Gaussians sitting on top of a polynomial continuum




Peaks can also be “data defined” meaning their amplitude is the area between the continuum observed counts

Recommend sticking with Gaussian + Polynomial unless bad fit



Peaks (continued)




Peaks can share a continuum

- If InterSpec doesn't automatically have peaks share a continuum: right-click on a peak, and choose to "Combine Cont. Right/Left" or "Own Continuum".

Peak continuum types

You can select different types of continuums by right-clicking on peak, or in “Peak Editor”



Peak Editor

Parameter	Value	Uncertainty	Fit
Centroid	383.4982	0.0037	<input checked="" type="checkbox"/>
FWHM	1.5922	0.0072	<input checked="" type="checkbox"/>
Amplitude	37988.8704	202.1275	<input checked="" type="checkbox"/>
χ^2/DOF	1.5639	0.0000	
ROI Start (keV)	379.4738	0.0000	
ROI End (keV)	387.5326	0.0000	
Cont. P0	546.2095	17.8410	<input checked="" type="checkbox"/>
Cont. P1	-37.9476	3.0323	<input checked="" type="checkbox"/>
Peak Color	0.0000	0.0000	
Nuclide			
Photopeak			
Label			
Peak Color	None		
Peak Type	Constant		
Continuum	<input checked="" type="checkbox"/> Linear		
Skew Type	Quadratic		
	Cubic		
	Flat Step		
	Linear Step		
	Bi-linear Step		
	Global Cont.		

Buttons: ? Cancel Accept Delete

More detailed and advanced peak-fitting will be covered in part two






Nuclide ID

general approach



Nuclide ID




The “Nuclide Search” tab lets you search on one or more energies

You can enter an energy and search window manually, or by clicking on a peak, or the spectrum

By default search-results ordered by “Profile” which generally is the best, but can also be ordered by any table column

Clicking on a result row in table will show reference lines for that nuclide, on the spectrum


Nuclide ID (continued)



Right-clicking on a peak after fitting it will also show you some suggestions

- Suggestions are meant to assist you – they are not an automated identification
- Especially useful for catching sum-peaks, or escape peaks

Nuclide ID (continued):



The "Reference Photopeaks" tab (covered next) is also really useful for nuclide ID



Reference Photopeaks

Some basics



Reference Photopeaks



The “Reference Photopeaks” tab might be the most useful tool in InterSpec

- Gives you access to all >4000 nuclides, all x-rays, and ~200 reactions
- You can age nuclides, add shielding, DRF, and when appropriate just prompt products
- The reference lines on the spectrum provide additional information when you mouse-over them

Reference Photopeaks (continued)




You can enter “nuclides” in fairly flexible, e.g., all the following are valid:

- Co56, Co-56, 56Co, 56 Co, 56-Co, cobalt 56, 56 Cobalt.
- Ho166m, Ho-166 meta, 166m Ho, 166m-Ho, holmium 166m
- Once you finish entering, will be converted to format like “Co56”

Reactions are entered like: “H(n,g)”, “Ge(n,n)”, etc

You can also enter “Background”, or an element name/symbol for its fluorescence x-rays


Reference Photopeaks (cont.)



The height of the lines on the spectrum indicate relative intensity of gammas/x-rays

- The heights are scaled so the highest intensity gamma line in currently displayed energy range has a height of 1, and all other lines are scaled linearly relative to this
- Even if the spectrum y-axis is log, reference lines are linear – this best matches what is seen in the data-spectrum, given the continuum counts
- Shielding and DRF (if entered) are taken into account in the chart-line amplitudes, but not on table values


Reference Photopeaks (cont.)



To display lines for multiple isotopes, click “Add Another” and enter a new nuclide in

- Each nuclide can have a different shielding, age, etc. Changing these fields only changes the current nuclide, if you want to alter a previously entered set of reference lines, just type that name into the “Nuclide” field
- Line height scaling for currently visible energy range happens on a per-nuclide level
- If you fit a peak while showing multiple nuclides, usually the best one is chosen to associate with a peak

Reference Photopeaks (cont.)




Reference lines for nuclides will also include x-rays that are emitted during decays

- x-ray intensities are computed during decay calculations so their intensities are comparable to the gammas
- But be aware of geometry or fluorescence effects


X-rays are shown as dotted lines on the spectrum

If you just enter a element, its fluorescence x-rays are shown

Reference Photopeaks (cont.)



Ir-192 reference lines with no shielding




Ir-192 reference lines with heavy shielding

Adding shielding can be useful to:


- Knock down high-intensity low-energy lines that are easily shielded
- Account for large amounts of shielding



Energy Calibration



Energy Calibration



The “Energy Calibration” tab provides extensive energy calibration capabilities

Fitting energy calibration coefficients, using peaks associated with nuclides is particularly useful

Will cover energy calibration more next week – you can also click (?) in lower-left of tab




Fitting activity and shielding based on peaks in the spectrum

A quick introduction - not covering theory, or many,
many details





Fitting activity, shielding, and age in InterSpec



First: fit peaks, associating them with nuclides

- Peaks that are in both the background and foreground do not need to be fit. Or if a peak is for a nuclide you don't care about you can skip it
- If you have lots of peaks for a nuclide, like a HPGe spectrum of Ra-226, you can maybe fit for just the larger ones

Fitting activity, shielding, and age in InterSpec (cont.)




From the "Tools" menu, select "Activity/Shielding Fit"

Activity/Shielding Fit Tool



Get an idea of how good the fit is



You can select peaks to use in the fit

Some options

Select DRF

Enter distance

Select shielding(s)


Fit for the quantities you want

Decide if you want activity and/or age fixed for each nuclide

Activity/Shielding Fit Tool (peaks to use)

All of the peaks you have fit for in the foreground will appear in this list.

- You can choose if a peak should participate in the fitting process.



Any nuclides associated with the peaks you have selected, will appear here

- You can manually edit an activity, and choose not to fit it, if you happen to know it, but not the shielding
- Not shown above, but age of nuclides is also editable, or fittable for applicable nuclides

Activity/Shielding Fit Tool (shielding)

Detector: LaBr 10%

Distance: 100 cm

Add Shielding: Material Generic

Fe (iron)

$\rho=7.9 \text{ g/cm}^3, \overline{AN}\sim 26.0$

Thickness 6.35 (± 0.35) mm Fit



Generic

AN 26 Fit AD 0.0 g/cm² Fit

Clicking on these icons will toggle between shielding types

- The shielding can be either a specific type of material from InterSpecs database, or it can be a generic material where you specify the atomic number (AN) and areal density (AD), or you can enter a chemical formula and density.
- You can choose to fit for, or fix, thickness, AN, and AD
- You can have multiple shielding's, but of course this can get degenerate to fit for multiple of them, especially if they have similar effective atomic numbers

Activity/Shielding Fit Tool: Graphic




The graphic in this tool is to help provide you an idea of how well the data lines up to the model

- Left: chart shows you how many statistical deviations each peak is away from the best-estimate answer
- Right: chart shows you the multiple of how far away from the model each peak-area is. The error bars are 1-sigma statistical uncertainties.
 - E.x., if the dot corresponding to a peak is at 2.0, that indicates that that peak would give double the activity that is currently showing

Generally, if after doing a "Perform Model Fit," a peak is greater than maybe 5-or-so-sigma away from the nominal value, you should consider removing that peak from the fit, or checking that its associated with the correct gamma line, or that there isn't an interfering nuclide that you aren't currently fitting for

Activity/Shielding Fit Tool: Detector Response Functions



If you click/tap on the detector area, it will bring up this dialog

- InterSpec includes a few generic Detector Response Functions (DRF)
- You can load DRFs from a few different common sources
- You can enter a mathematical formula for efficiency as a function of energy

The “Make Detector Response” tool can also be used to create your own DRF from characterization data

Generally the same DRF can be used for all detectors of the same model

Detector response functions not further covered in this presentation

Activity/Shielding Fit Tool – some words of caution



This tool really does require your critical thinking and spectroscopy knowledge.

You can shoot yourself in the foot.

And its also not magic. You can only fit for as many quantities as you have peaks for. So for Cs137, you can only fit for activity or shielding thickness; not both, or not the atomic number of shielding.

If you are fitting for the atomic number of shielding, you generally will need one or more peaks below 300 keV.


And you always need to ask yourself if the result makes sense.

And you should always add additional uncertainties for things like distance, or the detector response function, or a million other things.

Example Problem setup:

A 20cm x 20cm x 20cm package is found to be radioactive

A 20 minute spectrum is taken, 40 cm from surface of package using HPGe detector




Please: open "example_problem_1.n42", perform nuclide identification and quantification

Assume Lead (Pb) shielding, and source is centered in box

For convenience, the detector response function and background spectrum is included in the N42 file



Spectrum to be evaluated (only contains Eu152)



True Answer:

Eu-152

Result: 449 ± 6 kBq (12.1 ± 0.2 μ Ci)

6.3 ± 0.1 mm Pb shielding

Truth: 451 ± 3 kBq (12.2 ± 0.1 μ Ci)

7.0 mm Pb shielding

Uncertainties are $1\text{-}\sigma$, statistical only.

